Hemoflow F – Series

High Performance Steam

Capillary Dialysers

CE 0123

		F4 HPS	F5 HPS	F6 HPS	F7 HPS	F8 HPS	F10 HPS
Art.-No.		5007041	5007051	5007061	5007071	5007081	5007201
Q_B [mL/min]		200	200 / 300	200 / 300	200 / 300	300 / 400	300 / 400
Clearances (mL/min)							
$Q_D = 500$ mL/min, $UF = 0$	Urea	170	179 / 227	186 / 243	188 / 247	252 / 290	259 / 300
	Creatinine	149	162 / 196	173 / 215	175 / 220	224 / 251	230 / 259
	Phosphate	123	139 / 162	148 / 175	155 / 186	193 / 212	208 / 231
	Vitamin B ₁₂	75	84 / 91	92 / 100	102 / 113	118 / 124	131 / 139
K ₀ A Urea	mL/min	495	606	746	789	848	945
UF-coefficient	mL/h/mmHg	8	10	13	16	18	21
V (blood priming volume)	mL	51	63	78	96	113	132
Max. dialysate flow	mL/min	800	1000				
ΔP (pressure drop blood, $Q_B = 300$ mL/min)	mmHg	142	111	91	71	62	52
ΔP (pressure drop dialysate, $Q_D = 500$ mL/min)	mmHg	23	21	19	16	14	12
A (effective surface area)	m ²	0.8	1.0	1.3	1.6	1.8	2.2
TMP max.	mmHg	600					
Membrane		Fresenius Polysulfone®					
Sterilisation method		INLINE steam					

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
F-Series High Performance Steam Capillary Dialysers

GENERAL INFORMATION

Refer to product or carton label for:

	Single use only		Blood pathway sterile. Steam sterilised
	Expiry date		Refer to instructions for use
	Batch		Order number
	Date of manufacture		Storage temperature range
	Units		Use only on machines with exact UF-control

Indications: Hemoflow dialysers are designed for single use in chronic haemodialysis.

Contraindications: Contraindications are unknown. Generally, contraindications for haemodialysis are applicable.

Side-effects: In rare cases hypersensitivity reactions may occur during haemodialysis treatment. In severe cases dialysis must be discontinued and the appropriate medication initiated.

The dialyser is steam sterilised and thus contains no sterilisation residues. For patients with a known ethylene oxide hypersensitivity **all** components of the extracorporeal circuit should be ethylene oxide-free.

Anticoagulation: It is recommended to introduce an anticoagulant to the extracorporeal circuit. Nature, amount and method of application of an anticoagulant must be prescribed by the responsible physician (e.g. initial heparin bolus of 2000 IU followed by a continuous dose of 1000 IU/h up until the last hour of treatment). Coagulation should be monitored by a standard clotting time test.

Materials: Membrane: Fresenius Polysulfone®, Housing: Polycarbonate, Potting material: Polyurethane, O-ring: Silicone.

Further information may be obtained on request.

WARNINGS

Due to the relatively high hydraulic permeability of the membrane, the dialysers must only be used with dialysis machines which enable precise volumetric fluid control. Refer to dialysis machine manufacturer's technical and safety instructions for the use of high permeability membranes.

Use only if unit package is intact, sealing caps are in place and the dialyser is undamaged.

Dialysers must not be used after expiry date (see label).

Each dialyser is checked for integrity prior to leaving the factory. If a blood leak should arise, the dialyser must be exchanged.

The dialyser is intended for single-use only. Re-use may be hazardous to both the patient and operator. Cleansing solutions and disinfectants may damage materials employed for the housing, potting and membrane. Safety of use can no longer be guaranteed and the manufacturer assumes no liability.

HAEMODIALYSIS

Priming

Affix dialyser in a vertical position. Aseptically connect the arterial and venous blood lines to the lower and upper dialyser ports respectively. Attach the dialysate tubing connectors with the inlet to the upper and the outlet to the lower dialysate port. Fill dialyser according to the instructions provided with the machine (thereby the dialyser must be turned 180°). We recommend bag-free priming with ONLINE plus™ (priming volume: 500 mL). Alternatively the blood compartment can be filled and de-aerated in recirculation mode with isotonic saline solution (e.g. 500 mL bag). Ensure complete de-aeration of the blood compartment.

Patient Connection

Connect arterial line to the patient's circulation. Allow blood to flow in the blood line and dialyser (pump speed approx. 100 mL/min). Connect venous blood line to the patient.

Adjust treatment parameters to desired setting.

The weight loss should be monitored.

Termination

Perform reinfusion of the blood at the end of the treatment according to the instructions provided with the machine. We recommend the bag-free ONLINE plus™ procedure. Alternatively isotonic saline solution (e.g. 500 mL bag) can be used for reinfusion. The blood should be reinfused completely.

WARRANTY

Products with manufacturing defects will be replaced if defect is reported with details of the lot number.

The manufacturer will not be liable for any misuse, improper handling, non compliance with instructions for use and cautionary notes and for any damage incurred subsequent to the manufacturer's delivery of the dialyser.

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Fresenius Medical Care AG & Co. KGaA